

TRADITIONAL HIGH BACK CHAIRS WITH TORSO ARM

From **£250.00** Exc. VAT



HAN9120

HAN9130

Traditional High Back Chairs with Torso Arm at a glance:

- Pneumatic height adjustment – standard 20”-28”
- Pulled seat prevents upper leg impingement
- 360° swivel action
- Anti-microbial vinyl covering
- Hard free-wheeling castor wheels
- Contoured back support with dual adjustment
- Black base with hard free-wheeling castor wheels
- Height adjustable torso/arm
- Available with or without a footrest ring

Traditional High Back Chairs with Torso Arm dimensions (standard):

Low and High height versions are available.

- Height: 20” – 28” (Adjustable)
- Seat Pad Diameter: 15”
- Wheel Base: 23”
- Max. User Weight 165kg

Production lead time: approx. 7-10 working days

 UK Made. CE marked  Anti-microbial vinyl  Assembly required

 **FREE Delivery** anywhere in mainland UK by courier

WITHOUT FOOT REST

HAN9120 | **£250.00** Exc. VAT

WITH FOOT REST

HAN9130 | **£280.00** Exc. VAT

These Traditional High Back Chairs with Torso Arms feature comfortable, contoured back support with dual adjustment, a pneumatic height adjustment, a black base with hard, free-wheeling castor wheels and a full 360° range of motion. These luxury high-back chairs are available in a choice of versatile and hardwearing anti-microbial vinyl upholstery colours. Lastly, this traditional dentist's chair features a height-adjustable torso/arm that can provide comfort and support during time-intensive or intricate dental procedures.

This product is available with or without an integrated footrest ring.

Extra options:



Additional Memory Foam padding gives extra comfort if seated for long periods

+ £30



Wheel Variants:
Soft Anti-Static –roll freely on hard floors; Brake Loaded –roll & lock when sat on; Static –cup feet

+ £25



Height Variants:
Low Lift 19”-24”
High Lift 22”-32”

+ £20

Order online at meditelle-dental.co.uk or email sales@meditelleltd.co.uk

TRADITIONAL HIGH BACK CHAIRS WITH TORSO ARM

Specifications



HAN9120




HAN9130

UK Made. CE marked Anti-microbial vinyl Assembly required

FREE Delivery anywhere in mainland UK by courier

Technical information:



Height Information
 Standard 20"-28" (51cm - 71cm)
 Low option 19"-24" (48cm - 61cm)
 High option 22"-32" (56cm - 81cm)

SWL 165kg

Vinyl colours:



Special features:



Foot Rest Ring - provides a foot rest support for added comfort and if working at a higher height.

(HAN9130)



Torso Arm Rest - Height adjustable swing support.



Contoured moulded back support - features dual height and tilt adjustment with backrest pivot hinge to provide free float tilting for lumbar comfort.