

DENTAL OPERATOR CHAIRS

From **£180.00** Exc. VAT



HAN9240

HAN9220

Dental Operator Chairs at a glance:




- Pneumatic height adjustment – standard 20”-28”
- Contoured back support with dual adjustment
- 360° swivel action
- Anti-microbial vinyl covering
- Pulled seat prevents upper leg impingement
- Hard free-wheeling castor wheels
- Available with and without a footrest ring

Dental Operator Chair dimensions (standard):

Low and High height versions are available.

- Height: 20” – 28” (Adjustable)
- Seat Pad Diameter: 14”
- Wheel Base: 23”
- Max. User Weight 165kg

Production lead time: approx. 7-10 working days

 UK Made. CE marked  Anti-microbial vinyl  Assembly required

 **FREE Delivery** anywhere in mainland UK by courier

WITHOUT FOOT REST

HAN9240 | **£180.00** Exc. VAT

WITH FOOT REST

HAN9220 | **£210.00** Exc. VAT

Our Dental Operator Chairs with contoured, moulded back supports are the perfect addition to any busy dentist's office. Work in full comfort with a full height and tilt adjustment that provides maximum lumbar support. Meditelle's Dental Operator Chairs offer a smooth action pneumatic height adjustment (20” – 28”), a full 360° swivel range of motion and a choice of hardwearing anti-microbial vinyl colour options.

This practical Dental Operator Chair further benefits from hard, free-wheeling castor wheels and a pulled seat to help prevent upper leg impingement.

Our premium Dental Operator Chairs are available with or without an additional footrest ring.

Extra options:



Additional Memory Foam padding gives extra comfort if seated for long periods

+ £30



Wheel Variants:
Soft Anti-Static –roll freely on hard floors; Brake Loaded –roll & lock when sat on; Static –cup feet

+ £25



Height Variants:
Low Lift 19”-24”
High Lift 22”-32”

+ £20



Larger Seat Pad: increases overall seat pad diameter by +2”

+ £30

Order online at meditelle-dental.co.uk or email sales@meditelleltd.co.uk

DENTAL OPERATOR CHAIRS

Specifications



HAN9240




HAN9220

UK Made. CE marked Anti-microbial vinyl Assembly required

FREE Delivery anywhere in mainland UK by courier

Technical information:



Height Information
 Standard 20"-28" (51cm - 71cm)
 Low option 19"-24" (48cm - 61cm)
 High option 22"-32" (56cm - 81cm)

SWL 165kg

Vinyl colours:



Special features:



Foot Rest Ring - provides a foot rest support for added comfort and if working at a higher height.

(HAN9220)



Contoured moulded back support - features dual height and tilt adjustment with backrest pivot hinge to provide free float tilting for lumbar comfort.