

# DENTAL OPERATOR STOOLS WITH TORSO ARM

From **£220.00** Exc. VAT



HAN9270

HAN9280

## Dental Operator Stool with Torso Arm at a glance:

- Pneumatic height adjustment – standard 20”-28”
- Height adjustable torso/arm
- 360° swivel action
- Anti-microbial vinyl covering
- Pulled seat prevents upper leg impingement
- Hard free-wheeling castor wheels
- Available with and without a footrest ring


## Dental Operator Stool with Torso Arm dimensions (standard):

*Low and High height versions are available.*

- Height: 20” – 28” (Adjustable)
- Seat Pad Diameter: 14”
- Wheel Base: 23”
- Max. User Weight 165kg

**Production lead time: approx. 7-10 working days**

 UK Made. CE marked  Anti-microbial vinyl  Assembly required

 **FREE Delivery** anywhere in mainland UK by courier

**WITHOUT FOOT REST**

HAN9270 | **£220.00** Exc. VAT

**WITH FOOT REST**

HAN9280 | **£250.00** Exc. VAT

This Dental Operator Stool with Torso Arm is one of our most popular products. These elegant and high-quality dental chairs feature a smooth action pneumatic height adjustment (20” – 28”) and a practical 360° swivel for ease of motion. They also benefit from a height-adjustable torso/arm making them ideal for examinations and precision dental work.

Our Dental Operator Stool with Torso Arm are available in a variety of anti-microbial vinyl finishes and benefits from a pulled seat which helps to prevent upper leg impingement. Lastly, these versatile dentist’s chairs benefit from an easy clean base and hard, free-wheeling castor wheels for manoeuvrability.

*This popular swivel arm dentist’s stool is available with or without an additional footrest ring for added leg support!*

## Extra options:



Additional Memory Foam padding gives extra comfort if seated for long periods

**+ £30**



Wheel Variants:  
Soft Anti-Static –roll freely on hard floors; Brake Loaded –roll & lock when sat on; Static –cup feet

**+ £25**



Height Variants:  
Low Lift 19”-24”  
High Lift 22”-32”

**+ £20**



Larger Seat Pad: increases overall seat pad diameter by +2”

**+ £30**

Order online at [meditelle-dental.co.uk](https://meditelle-dental.co.uk) or email [sales@meditelleltd.co.uk](mailto:sales@meditelleltd.co.uk)

# DENTAL OPERATOR STOOLS WITH TORSO ARM

## Specifications




HAN9270

HAN9280

UK Made. CE marked Anti-microbial vinyl Assembly required

**FREE Delivery** anywhere in mainland UK by courier

### Technical information:



**Height Information**  
 Standard 20"-28" (51cm - 71cm)  
 Low option 19"-24" (48cm - 61cm)  
 High option 22"-32" (56cm - 81cm)

SWL 165kg

### Vinyl colours:



### Special features:



**Foot Rest Ring** - provides a foot rest support for added comfort and if working at a higher height.

(HAN9280)



**Torso Arm Rest** - Height adjustable swing support.

Order online at [meditelle-dental.co.uk](http://meditelle-dental.co.uk) or email [sales@meditelleltd.co.uk](mailto:sales@meditelleltd.co.uk)