

TRI-SHAPED STOOLS WITH TORSO ARM

From **£295.00** Exc. VAT



HAN9330

HAN9335

Tri-Shaped Stools with Torso Arm at a glance:

- Pneumatic height adjustment – standard 20”-28”
- Tri-shaped seat prevents upper leg impingement
- Height adjustable torso/arm
- 360° swivel action
- Anti-microbial vinyl covering
- Hard free-wheeling castor wheels
- Available with and without a footrest ring

Tri-Shaped Stool with Torso Arm dimensions (standard):

Low and High height versions are available.

- Height: 20” – 28” (Adjustable)
- Seat Pad Diameter: 15”
- Wheel Base: 23”
- Max. User Weight 165kg

Production lead time: approx. 7-10 working days



UK Made. CE marked



Anti-microbial vinyl



Assembly required



FREE Delivery anywhere in mainland UK by courier

WITHOUT FOOT REST

HAN9330 | **£295.00** Exc. VAT

WITH FOOT REST

HAN9335 | **£325.00** Exc. VAT

Our Tri Shaped Stools with Torso Arms are built with comfort and stability in mind making them the ideal seating solution for any dentist's office or surgery. These functional tri-shaped stools boast a tri-shaped seat pad as standard, which is shaped specifically to prevent upper leg impingement while hard, free-wheeling castors keep you firmly rooted to the floor.

This tri-shaped stool with a torso arm is available in a wide range of antimicrobial vinyl colour options and boasts a pneumatic height adjustment and 360° swivel action as well as an adjustable torso/arm which is ideally suited for dentists and surgeons needing to perform precision or time-intensive procedures.

Our Meditelle Tri-Shaped Stools with Torso Arms are available with and without a footrest ring.

Extra options:



Additional Memory Foam padding gives extra comfort if seated for long periods

+ £30



Wheel Variants:
Soft Anti-Static -roll freely on hard floors; Brake Loaded -roll & lock when sat on; Static -cup feet

+ £25



Height Variants:
Low Lift 19”-24”
High Lift 22”-32”

+ £20

Order online at meditelle-dental.co.uk or email sales@meditelleltd.co.uk

TRI-SHAPED STOOLS WITH TORSO ARM

Specifications



HAN9330



HAN9335



UK Made. CE marked



Anti-microbial vinyl




Assembly required



FREE Delivery anywhere in mainland UK by courier

Technical information:



Height Information
Standard 20"-28" (51cm - 71cm)
Low option 19"-24" (48cm - 61cm)
High option 22"-32" (56cm - 81cm)



SWL 165kg

Vinyl colours:

	Burnt Red		Rose		Mulberry		Lilac		Navy
	Royal		Lido		Aqua		Ocean		Citrus
	Sage		Grey		Dove		Ginger		Buttercup
	Taupe		Barley		Stone		White		Black

Special features:



Foot Rest Ring - provides a foot rest support for added comfort and if working at a higher height.

(HAN9335)



Torso Arm Rest - Height adjustable swing support.