

TRI-SHAPED CHAIRS WITH TORSO ARM

From **£335.00** Exc. VAT



HAN9340

HAN9345

Tri-Shaped Chairs with Torso Arm at a glance:

- Pneumatic height adjustment – standard 20”-28”
- Height adjustable torso/arm
- Tri-shaped seat prevents upper leg impingement
- 360° swivel action
- Anti-microbial vinyl covering
- Hard free-wheeling castor wheels
- Contoured back support with dual adjustment
- Available with and without a footrest ring

Tri-Shaped Chair with Torso Arm dimensions (standard):

Low and High height versions are available.

- Height: 20” – 28” (Adjustable)
- Seat Pad Diameter: 15”
- Wheel Base: 23”
- Max. User Weight 165kg

Production lead time: approx. 7-10 working days

 UK Made. CE marked  Anti-microbial vinyl  Assembly required

 **FREE Delivery** anywhere in mainland UK by courier

WITHOUT FOOT REST

HAN9340 | **£335.00** Exc. VAT

WITH FOOT REST

HAN9345 | **£365.00** Exc. VAT

Our Tri-Shaped Chairs with Torso Arms are built with comfort and stability in mind making them the ideal seating solution for any dentist's office or surgery. These functional tri-shaped chairs boast a tri-shaped seat pad as standard, which is shaped specifically to prevent upper leg impingement, hard, free-wheeling castors keep you firmly rooted to the floor and a contoured back support with dual adjustment for maximum comfort and posture management.

This tri-shaped chair is available in a wide range of antimicrobial vinyl colour options and boasts a pneumatic height adjustment and 360° swivel action as well as a practical, height-adjustable torso/arm for added support.

Our Meditelle Tri Shaped Chairs with Torso Arms are available with and without a footrest ring.

Extra options:



Additional Memory Foam padding gives extra comfort if seated for long periods

+ £30



Wheel Variants:
Soft Anti-Static –roll freely on hard floors; Brake Loaded –roll & lock when sat on; Static –cup feet

+ £25



Height Variants:
Low Lift 19”-24”
High Lift 22”-32”

+ £20

Order online at meditelle-dental.co.uk or email sales@meditelleltd.co.uk

TRI-SHAPED CHAIRS WITH TORSO ARM

Specifications




HAN9340

HAN9345

UK Made. CE marked Anti-microbial vinyl Assembly required

FREE Delivery anywhere in mainland UK by courier

Technical information:



Height Information
 Standard 20"-28" (51cm - 71cm)
 Low option 19"-24" (48cm - 61cm)
 High option 22"-32" (56cm - 81cm)

SWL 165kg

Vinyl colours:



Special features:



Foot Rest Ring - provides a foot rest support for added comfort and if working at a higher height.

(HAN9345)



Torso Arm Rest - Height adjustable swing support.