

# TILT SADDLE STOOLS WITH TORSO ARM

From **£300.00** Exc. VAT



HAN9440

HAN9445

## Tilt Saddle Stools with Torso Arm at a glance:

- Pneumatic height adjustment – standard 20”-28”
- Contoured saddle-shaped pad aids in correct posture
- Pneumatic 35° seat tilt adjustment
- Height adjustable torso/arm with 270° rotation
- 360° swivel action
- Anti-microbial vinyl or black leather covering
- Hard free-wheeling castor wheels
- Available with and without a footrest ring

## Tilt Saddle Stool with Torso Arm dimensions (standard):

*Low and High height versions are available.*

- Height: 20” – 28” (Adjustable)
- Wheel Base: 25”
- Max. User Weight 165kg

**Production lead time: approx. 7-10 working days**

- UK Made. CE marked
- Anti-microbial vinyl
- Assembly required
- FREE Delivery** anywhere in mainland UK by courier

**WITHOUT FOOT REST**

HAN9440 | **£300.00** Exc. VAT

**WITH FOOT REST**

HAN9445 | **£330.00** Exc. VAT

Introducing our maximum comfort tilt saddle stools with torso arm, these luxury stools adopt a new and ergonomic style of seating in the form of a contoured saddle-shaped seat pad which provides expert comfort and helps to correct posture when used long-term, in addition, this version of the stool benefits from a height adjustable, rotating torso arm with 270° rotation which can be used for aiding precision dental work.

These dentist's saddle stools with torso arm offer a pneumatic height adjustment and a 360° swivel action for maximum manoeuvrability. This saddle stool is available in a wide range of anti-microbial vinyl colour options with an added option of a black leather covering and boasts hard free-wheeling castor wheels for added stability.

Our tilt saddle stools with torso arm are available with and without an additional footrest ring which provides added comfort when working at height. Choose the option that's right for you.

## Extra options:



Additional Memory Foam padding gives extra comfort if seated for long periods

**+ £30**



Wheel Variants:  
Soft Anti-Static –roll freely on hard floors; Brake Loaded –roll & lock when sat on; Static –cup feet

**+ £25**




Height Variants:  
Low Lift 19”-24”  
High Lift 22”-32”

**+ £20**



Smaller Saddle Pad: reduces overall seat pad diameter by -1”

**+ £30**



Alternative Black Leather Upholstery

**+ £65**

Order online at [meditelle-dental.co.uk](https://meditelle-dental.co.uk) or email [sales@meditelleltd.co.uk](mailto:sales@meditelleltd.co.uk)

# TILT SADDLE STOOLS WITH TORSO ARM

## Specifications




HAN9440

HAN9445

UK Made. CE marked
 Anti-microbial vinyl
 Assembly required

**FREE Delivery** anywhere in mainland UK by courier

### Technical information:



### Vinyl colours:



### Special features:



**Foot Rest Ring** - provides a foot rest support for added comfort and if working at a higher height.

(HAN9445)



**Tilting Seat Pad** with 35° pneumatic seat tilt adjustment.



**Torso Arm Rest** - Height adjustable swing support.