

# TILT SADDLE CHAIR WITH TORSO ARM

From **£385.00** Exc. VAT



HAN9460

HAN9465

### Tilt Saddle Chairs with Torso Arm at a glance:

- Pneumatic height adjustment – standard 20”-28”
- A contoured saddle-shaped pad aids in correct posture
- Pneumatic 35° seat tilt adjustment
- 360° swivel action
- Anti-microbial vinyl covering
- Hard free-wheeling castor wheels
- Contoured back support with dual adjustment
- Height adjustable torso/arm with 270° rotation
- Available with and without a footrest ring

### Tilt Saddle Chair with Torso Arm dimensions (standard):

Low and High height versions are available.

- Height: 20” – 28” (Adjustable)
- Wheel Base: 25”
- Max. User Weight 165kg

**Production lead time: approx. 7-10 working days**

UK Made. CE marked Anti-microbial vinyl Assembly required

**FREE Delivery** anywhere in mainland UK by courier

**WITHOUT FOOT REST**

HAN9460 | **£385.00** Exc. VAT

**WITH FOOT REST**

HAN9465 | **£415.00** Exc. VAT

Our Tilt Saddle Chairs with Torso Arms are the perfect seating option for dentists who need to perform regular precision work. This premium dental chair features a contoured saddle-shaped seating pad which helps to aid correct posture, a 360° range of motion and a pneumatic 35° seat tilt adjustment. On top of this, this chair benefits from a height-adjustable torso/arm with 270° rotation, pneumatic height adjustment, contoured back support with dual adjustment and hard, free-wheeling castors for added stability.

Our Tilt Saddle Chairs with Torso Arm are available in a range of anti-microbial vinyl colours and an optional foot-rest ring which provides increased comfort when working with the seat fully extended.

### Extra options:



Additional Memory Foam padding gives extra comfort if seated for long periods

**+ £30**



Wheel Variants:  
Soft Anti-Static –roll freely on hard floors; Brake Loaded –roll & lock when sat on; Static –cup feet

**+ £25**




Height Variants:  
Low Lift 19”-24”  
High Lift 22”-32”

**+ £20**



Smaller Saddle Pad: reduces overall seat pad diameter by -1”

**+ £30**



Alternative Black Leather Upholstery

**+ £65**

Order online at [meditelle-dental.co.uk](https://meditelle-dental.co.uk) or email [sales@meditelleltd.co.uk](mailto:sales@meditelleltd.co.uk)

# TILT SADDLE STOOLS WITH TORSO ARM

## Specifications




HAN9460

HAN9465

UK Made. CE marked Anti-microbial vinyl Assembly required

**FREE Delivery** anywhere in mainland UK by courier

### Technical information:



### Vinyl colours:



### Special features:



**Foot Rest Ring** - provides a foot rest support for added comfort and if working at a higher height.

(HAN9465)



**Tilting Seat Pad** with 35° pneumatic seat tilt adjustment.



**Torso Arm Rest** - Height adjustable swing support.



**Contoured moulded back support** - features dual height and tilt adjustment with backrest pivot hinge to provide free float tilting for lumbar comfort.