

SURGEONS FOOT OPERATED TILT SADDLE CHAIRS

From **£330.00** Exc. VAT



HAN9485

HAN9475

Surgeons Foot Operated Tilt Saddle Chairs at a glance:

- Pneumatic foot operated height adjustment – standard 22"–27"
- A contoured saddle-shaped pad aids in correct posture
- Contoured back support with dual adjustment
- Pneumatic 35° seat tilt adjustment
- 360° swivel action
- Anti-microbial vinyl covering or black leather covering
- Hard free-wheeling castor wheels
- Available with height adjustable torso/arm with 270° rotation


Surgeons Foot Operated Tilt Saddle Chair dimensions (standard):

Low and High height versions are available.

- Height: 22" – 27" (Adjustable)
- Wheel Base: 25"
- Max. User Weight 165kg

Production lead time: approx. 7-10 working days

 UK Made. CE marked  Anti-microbial vinyl  Assembly required

 **FREE Delivery** anywhere in mainland UK by courier

WITHOUT TORSO/ARM SUPPORT

HAN9485 | **£415.00** Exc. VAT

WITH TORSO/ARM SUPPORT

HAN9475 | **£515.00** Exc. VAT

Our Surgeon's Foot Operated Tilt Saddle Chairs feature a comfortable and ergonomic contoured saddle-shaped seat pad which helps to correct poor posture and create a healthy seating alternative to standard dentist's chairs. This premium saddle chair boasts a foot-operated pneumatic height adjustment, contoured back support with dual adjustment and a 360° swivel for maximum manoeuvrability. On top of this, this saddle chair benefits from a pneumatic 35° seat tilt adjustment, a choice of anti-microbial vinyl colour options or black leather covering and hard free-wheeling castors for stability.

This Surgeons Foot Operated Tilt Saddle Chair is available with and without an additional torso arm with 270° rotation. Choose the right option for you from the dropdown menu. The addition of an armrest helps to provide added stability and comfort when performing precision dental work.

Choose the perfect chair for your dental surgery at Meditelle Dental!

Extra options:



Additional Memory Foam padding gives extra comfort if seated for long periods

+ £30



Wheel Variants:
Soft Anti-Static –roll freely on hard floors; Brake Loaded –roll & lock when sat on; Static –cup feet

+ £25




Height Variants:
Low Lift 18"–22"
High Lift 23"–31"

+ £20



Smaller Saddle Pad:
reduces overall seat pad diameter by ~1"

+ £30



Alternative Black Leather Upholstery

+ £65

Order online at meditelle-dental.co.uk or email sales@meditelleltd.co.uk

SURGEONS FOOT OPERATED TILT SADDLE CHAIRS

Specifications



HAN9485



HAN9475



UK Made. CE marked



Anti-microbial vinyl




Assembly required



FREE Delivery anywhere in mainland UK by courier

Technical information:



Vinyl colours:

	Burnt Red		Rose		Mulberry		Lilac		Navy
	Royal		Lido		Aqua		Ocean		Citrus
	Sage		Grey		Dove		Ginger		Buttercup
	Taupe		Barley		Stone		White		Black

Special features:



Pneumatic Foot Operated Height Control – allows hands free, sterile height adjustment of the seat by foot.



Tilting Seat Pad with 35° pneumatic seat tilt adjustment.



Torso Arm Rest – Height adjustable swing support. (HAN9475)



Contoured moulded back support – features dual height and tilt adjustment with backrest pivot hinge to provide free float tilting for lumbar comfort.