

SURGEONS FOOT OPERATED TILT SADDLE STOOLS



SADDLE SEATS

From **£330.00** Exc. VAT



HAN9490

HAN9495

Surgeons Foot Operated Tilt Saddle Stools at a glance:

- Pneumatic foot operated height adjustment – standard 22"–27"
- A contoured saddle-shaped pad aids in correct posture
- Pneumatic 35° seat tilt adjustment
- 360° swivel action
- Anti-microbial vinyl or black leather covering
- Hard free-wheeling castor wheels
- Available with height adjustable torso/arm rest with 270° rotation

Surgeons Foot Operated Tilt Saddle Stool dimensions (standard):

Low and High height versions are available.

- Height: 22" – 27" (Adjustable)
- Wheel Base: 25"
- Max. User Weight 165kg

Production lead time: approx. 7-10 working days



UK Made. CE marked



Anti-microbial vinyl



Assembly required



FREE Delivery anywhere in mainland UK by courier

WITHOUT TORSO/ARM SUPPORT

HAN9490 | **£330.00** Exc. VAT

WITH TORSO/ARM SUPPORT

HAN9495 | **£430.00** Exc. VAT

Our Surgeon's Foot Operated Tilt Saddle Stools feature a comfortable and ergonomic contoured saddle-shaped seat pad which helps to correct poor posture and create a healthy seating alternative to standard dentist's stools. This premium saddle stool boasts a foot-operated pneumatic height adjustment and a 360° swivel for maximum manoeuvrability. On top of this, this saddle stool benefits from a pneumatic 35° seat tilt adjustment, a choice of anti-microbial vinyl colour options or black leather covering and hard free-wheeling castors for stability.

This Surgeons Foot Operated Tilt Saddle Stool is available with and without an additional torso/arm rest with 270° rotation. Choose the right option for you. The addition of an armrest helps to provide added stability and comfort when performing precision dental work.

Extra options:



Additional Memory Foam padding gives extra comfort if seated for long periods

+ £30



Wheel Variants:
Soft Anti-Static –roll freely on hard floors; Brake Loaded –roll & lock when sat on; Static –cup feet

+ £25



Height Variants:
Low Lift 18"–22"
High Lift 23"–31"

+ £20



Smaller Saddle Pad:
reduces overall seat pad diameter by ~1"

+ £30



Alternative Black Leather Upholstery

+ £65

Order online at meditelle-dental.co.uk or email sales@meditelleltd.co.uk

SURGEONS FOOT OPERATED TILT SADDLE STOOLS

Specifications




HAN9490

HAN9495

UK Made. CE marked Anti-microbial vinyl Assembly required

FREE Delivery anywhere in mainland UK by courier

Technical information:



Vinyl colours:

Burnt Red	Rose	Mulberry	Lilac	Navy
Royal	Lido	Aqua	Ocean	Citrus
Sage	Grey	Dove	Ginger	Buttercup
Taupe	Barley	Stone	White	Black

Special features:



Pneumatic Foot Operated Height Control – allows hands free, sterile height adjustment of the seat by foot.



Tilting Seat Pad with 35° pneumatic seat tilt adjustment.



Torso Arm Rest – Height adjustable swing support.
(HAN9495)