

MEDI SADDLE STOOLS



MEDI SADDLE SEATS

From **£250.00** Exc. VAT

Saddle seating trial available. Visit meditelle-dental.co.uk/saddle-trial



HAN9510

HAN9515

Medi Saddle Stools at a glance:

- Pneumatic height adjustment – standard 20"–28"
- Contoured saddle-shaped pad aids in correct posture
- Pneumatic 35° seat tilt adjustment
- 360° swivel action
- Anti-microbial vinyl covering or black leather covering
- Hard free-wheeling castor wheels
- Available with and without a footrest ring

Medi Saddle Stool dimensions (standard):

Low and High height versions are available.

- Height: 20" – 28" (Adjustable)
- Wheel Base: 25"
- Max. User Weight 165kg

Production lead time: approx. 10-15 working days



UK Made. CE marked



Anti-microbial vinyl



Assembly required



FREE Delivery anywhere in mainland UK by courier

WITHOUT FOOT REST

HAN9510 | **£250.00** Exc. VAT

WITH FOOT REST

HAN9515 | **£280.00** Exc. VAT

Introducing our maximum comfort medi saddle stools, these luxury stools adopt a new and ergonomic style of seating in the form of a contoured saddle-shaped seat pad which provides expert comfort and helps to correct posture when used long-term, the seating pad of this stool can also be fine-tuned thanks to a pneumatic 35° seat tilt adjustment.

These medi saddle stools offer a pneumatic height adjustment and a 360° swivel action for maximum manoeuvrability.

As with all the chairs and stools in our range, this medi saddle stool is available in various anti-microbial vinyl colour options with an added option of a black leather covering and boasts hard free-wheeling castor wheels for added stability.

Our medi saddle stools are available with and without an additional footrest ring, providing added comfort when working at height. Choose the option that's right for you.

Extra options:



Additional Memory Foam padding gives extra comfort if seated for long periods

+ £30




Wheel Variants:
Soft Anti-Static -roll freely on hard floors; Brake Loaded -roll & lock when sat on; Static -cup feet

+ £25



Height Variants:
Low Lift 19"–24"
High Lift 22"–32"

+ £20



Alternative Black Leather Upholstery

+ £65

Order online at meditelle-dental.co.uk or email sales@meditelleltd.co.uk

MEDI SADDLE STOOLS

Specifications



HAN9510




HAN9515

UK Made. CE marked
 Anti-microbial vinyl
 Assembly required

FREE Delivery anywhere in mainland UK by courier

Technical information:



Vinyl colours:

Burnt Red	Rose	Mulberry	Lilac	Navy
Royal	Lido	Aqua	Ocean	Citrus
Sage	Grey	Dove	Ginger	Buttercup
Taupe	Barley	Stone	White	Black

Special features:



Foot Rest Ring - provides a foot rest support for added comfort and if working at a higher height.

(HAN9515)



Tilting Seat Pad with 35° pneumatic seat tilt adjustment.