

# MEDI SADDLE STOOLS WITH TORSO ARM



MEDI SADDLE SEATS

From **£350.00** Exc. VAT

Saddle seating trial available. Visit [meditelle-dental.co.uk/saddle-trial](https://meditelle-dental.co.uk/saddle-trial)



HAN9540

HAN9545

## Medi Saddle Stools with Torso Arm at a glance:

- Pneumatic height adjustment – standard 20"–28"
- Contoured saddle-shaped pad aids in correct posture
- Pneumatic 35° seat tilt adjustment
- Height adjustable torso/arm with 360° rotation
- 360° swivel action
- Anti-microbial vinyl or black leather covering
- Hard free-wheeling castor wheels
- Available with and without a footrest ring

## Medi Saddle Stool with Torso Arm dimensions (standard):

*Low and High height versions are available.*

- Height: 20" – 28" (Adjustable)
- Wheel Base: 25"
- Max. User Weight 165kg

**Production lead time: approx. 10-15 working days**



UK Made. CE marked



Anti-microbial vinyl



Assembly required



**FREE Delivery** anywhere in mainland UK by courier

**WITHOUT FOOT REST**

HAN9540 | **£350.00** Exc. VAT

**WITH FOOT REST**

HAN9545 | **£380.00** Exc. VAT

Introducing our full comfort medi saddle stools with torso arm and pneumatic seat adjustment, these luxury stools adopt a new and ergonomic style of seating in the form of a contoured saddle-shaped seat pad which provides expert comfort and helps to correct posture when used long-term. This version of the medi stool benefits from a height adjustable, 360° rotating torso arm which can be used for aiding precision dental work.

These medi saddle stools with torso arm offer a pneumatic height adjustment, 35° seat tilt adjustment and a 360° swivel action for maximum manoeuvrability. As with all the chairs and stools in our range, this saddle stool with torso arm is available in various anti-microbial vinyl colour options with an added option of a black leather covering and boasts hard free-wheeling castor wheels for added stability.

Our medi saddle stools with torso arm are available with and without an additional footrest ring which provides added comfort when working at height. Choose the option that's right for you.

## Extra options:



Additional Memory Foam padding gives extra comfort if seated for long periods

**+ £30**




Wheel Variants:  
Soft Anti-Static -roll freely on hard floors; Brake Loaded -roll & lock when sat on; Static -cup feet

**+ £25**



Height Variants:  
Low Lift 19"–24"  
High Lift 22"–32"

**+ £20**



Alternative Black Leather Upholstery

**+ £65**

Order online at [meditelle-dental.co.uk](https://meditelle-dental.co.uk) or email [sales@meditelleltd.co.uk](mailto:sales@meditelleltd.co.uk)

# MEDI SADDLE STOOLS WITH TORSO ARM

## Specifications




HAN9540

HAN9545

UK Made. CE marked Anti-microbial vinyl Assembly required

**FREE Delivery** anywhere in mainland UK by courier

### Technical information:



### Vinyl colours:

Burnt Red	Rose	Mulberry	Lilac	Navy
Royal	Lido	Aqua	Ocean	Citrus
Sage	Grey	Dove	Ginger	Buttercup
Taupe	Barley	Stone	White	Black

### Special features:



**Foot Rest Ring** - provides a foot rest support for added comfort and if working at a higher height.

(HAN9545)



**Tilting Seat Pad** with 35° pneumatic seat tilt adjustment.



**Torso Arm Rest** - Height adjustable swing support.