

# MEDI SADDLE CHAIRS

From **£335.00** Exc. VAT



HAN9550



HAN9555

## Medi Saddle Chairs at a glance:




- Pneumatic height adjustment – standard 20”-28”
- A contoured saddle-shaped pad aids in correct posture
- Pneumatic 35° seat tilt adjustment
- 360° swivel action
- Anti-microbial vinyl covering
- Hard free-wheeling castor wheels
- Contoured back support with dual adjustment
- Available with and without a footrest ring

## Medi Saddle Chair dimensions (standard):

Low and High height versions are available.

- Height: 20” – 28” (Adjustable)
- Wheel Base: 25”
- Max. User Weight 165kg

**Production lead time: approx. 10-15 working days**

 UK Made. CE marked  Anti-microbial vinyl  Assembly required

 **FREE Delivery** anywhere in mainland UK by courier

**WITHOUT FOOT REST**

HAN9550 | **£335.00** Exc. VAT

**WITH FOOT REST**

HAN9555 | **£365.00** Exc. VAT

Our Medi Saddle Chairs are the perfect seating option for dentists who need to perform regular precision work. This premium dental chair features a contoured saddle-shaped seating pad which helps to aid correct posture, a 360° range of motion and a pneumatic 35° seat tilt adjustment. On top of this, this chair benefits from a pneumatic height adjustment, contoured back support with dual adjustment and hard, free-wheeling castors for added stability.

Our Medi Saddle Chairs are available in a range of anti-microbial vinyl colours and an optional foot-rest ring which provides increased comfort when working with the seat fully extended. An alternative leather upholstery is also available for this item!

## Extra options:



Additional Memory Foam padding gives extra comfort if seated for long periods

**+ £30**




Wheel Variants:  
Soft Anti-Static –roll freely on hard floors; Brake Loaded –roll & lock when sat on; Static –cup feet

**+ £25**



Height Variants:  
Low Lift 19”-24”  
High Lift 22”-32”

**+ £20**



Alternative Black Leather Upholstery

**+ £65**

Order online at [meditelle-dental.co.uk](https://meditelle-dental.co.uk) or email [sales@meditelleltd.co.uk](mailto:sales@meditelleltd.co.uk)

# MEDI SADDLE CHAIRS

## Specifications



HAN9550




HAN9555

UK Made. CE marked   Anti-microbial vinyl   Assembly required

**FREE Delivery** anywhere in mainland UK by courier

### Technical information:



### Vinyl colours:



### Special features:



**Foot Rest Ring** - provides a foot rest support for added comfort and if working at a higher height.

(HAN9555)



**Tilting Seat Pad** with 35° pneumatic seat tilt adjustment.



**Contoured moulded back support** - features dual height and tilt adjustment with backrest pivot hinge to provide free float tilting for lumbar comfort.

Order online at [meditelle-dental.co.uk](http://meditelle-dental.co.uk) or email [sales@meditelleltd.co.uk](mailto:sales@meditelleltd.co.uk)