

MEDI SADDLE CHAIRS WITH TORSO ARM



MEDI SADDLE SEATS

From **£435.00** Exc. VAT

Saddle seating trial available. Visit meditelle-dental.co.uk/saddle-trial



HAN9560

HAN9565

Medi Saddle Chairs with Torso Arm at a glance:

- Pneumatic height adjustment – standard 20"–28"
- A contoured saddle-shaped pad aids in correct posture
- Pneumatic 35° seat tilt adjustment
- 360° swivel action
- Anti-microbial vinyl or black leather covering
- Hard free-wheeling castor wheels
- Contoured back support with dual adjustment
- Height adjustable torso/arm with 360° rotation
- Available with and without a footrest ring

Medi Saddle Chair with Torso Arm dimensions (standard):

Low and High height versions are available.

- Height: 20" – 28" (Adjustable)
- Wheel Base: 25"
- Max. User Weight 165kg

Production lead time: approx. 10-15 working days



UK Made. CE marked



Anti-microbial vinyl



Assembly required



FREE Delivery anywhere in mainland UK by courier

WITHOUT FOOT REST

HAN9560 | **£435.00** Exc. VAT

WITH FOOT REST

HAN9565 | **£465.00** Exc. VAT

Our Medi Saddle Chairs with Torso Arms are the perfect seating option for dentists who need to perform regular precision work. This premium dental chair features a contoured saddle-shaped seating pad which helps to aid correct posture, a 360° range of motion and a pneumatic 35° seat tilt adjustment. On top of this, this chair benefits from a height-adjustable torso/arm with 360° rotation, pneumatic height adjustment, contoured back support with dual adjustment and hard, free-wheeling castors for added stability.

Our Medi Saddle Chairs with Torso Arm are available in a range of anti-microbial vinyl colours and an optional foot-rest ring which provides increased comfort when working with the seat fully extended. An alternative leather upholstery is also available for this item!

Extra options:



Additional Memory Foam padding gives extra comfort if seated for long periods

+ £30




Wheel Variants:
Soft Anti-Static -roll freely on hard floors; Brake Loaded -roll & lock when sat on; Static -cup feet

+ £25



Height Variants:
Low Lift 19"–24"
High Lift 22"–32"

+ £20



Alternative Black Leather Upholstery

+ £65

Order online at meditelle-dental.co.uk or email sales@meditelleltd.co.uk

MEDI SADDLE CHAIRS WITH TORSO ARM

Specifications



HAN9560



HAN9565



UK Made. CE marked



Anti-microbial vinyl




Assembly required



FREE Delivery anywhere in mainland UK by courier

Technical information:



Vinyl colours:



Special features:



Foot Rest Ring - provides a foot rest support for added comfort and if working at a higher height.

(HAN9565)



Tilting Seat Pad with 35° pneumatic seat tilt adjustment.



Torso Arm Rest - Height adjustable swing support.



Contoured moulded back support - features dual height and tilt adjustment with backrest pivot hinge to provide free float tilting for lumbar comfort.