

SURGEONS FOOT OPERATED MEDI SADDLE CHAIRS



MEDI SADDLE SEATS

From **£465.00** Exc. VAT



HAN9570

HAN9590

Surgeons Foot Operated Medi Saddle Chairs at a glance:

- Pneumatic foot operated height adjustment – standard 22”-27”
- A contoured saddle-shaped pad aids in correct posture
- Contoured back support with dual adjustment
- Pneumatic 35° seat tilt adjustment
- 360° swivel action
- Anti-microbial vinyl or black leather covering
- Hard free-wheeling castor wheels
- Available with height adjustable torso/arm with 360° rotation

Surgeons Foot Operated Medi Saddle Chair dimensions (standard):

Low and High height versions are available.

- Height: 22” – 27” (Adjustable)
- Wheel Base: 25”
- Max. User Weight 165kg

Production lead time: approx. 10-15 working days

UK Made. CE marked Anti-microbial vinyl Assembly required

FREE Delivery anywhere in mainland UK by courier

WITHOUT TORSO/ARM SUPPORT HAN9570 | **£465.00** Exc. VAT

WITH TORSO/ARM SUPPORT HAN9590 | **£565.00** Exc. VAT

Our Surgeon's Foot Operated Medi Saddle Chairs feature a comfortable and ergonomic contoured saddle-shaped seat pad which helps to correct poor posture and create a healthy seating alternative to standard dentist's chairs. This premium saddle chair boasts a foot-operated pneumatic height adjustment, contoured back support with dual adjustment and a 360° swivel for maximum manoeuvrability. On top of this, this saddle chair benefits from a pneumatic 35° seat tilt adjustment, a choice of anti-microbial vinyl colour options or black leather covering and hard free-wheeling castors for stability.

This Surgeons Foot Operated Tilt Saddle Chair is available with and without an additional torso arm with 360° rotation. Choose the right option for you. The addition of an armrest helps to provide added stability and comfort when performing precision dental work.

Get the perfect dental stool today at Meditelle Dental.

Extra options:



Additional Memory Foam padding gives extra comfort if seated for long periods

+ £30




Wheel Variants:
Soft Anti-Static -roll freely on hard floors; Brake Loaded -roll & lock when sat on; Static -cup feet

+ £25



Height Variants:
Low Lift 18”-22”
High Lift 23”-31”

+ £20



Alternative Black Leather Upholstery

+ £65

Order online at meditelle-dental.co.uk or email sales@meditelleltd.co.uk

SURGEONS FOOT OPERATED MEDI SADDLE CHAIRS

Specifications



HAN9570




HAN9590

UK Made. CE marked
 Anti-microbial vinyl
 Assembly required

FREE Delivery anywhere in mainland UK by courier

Technical information:



Vinyl colours:



Special features:



Pneumatic Foot Operated Height Control - allows hands free, sterile height adjustment of the seat by foot.



Tilting Seat Pad with 35° pneumatic seat tilt adjustment.



Torso Arm Rest - Height adjustable swing support. **(HAN9590)**



Contoured moulded back support - features dual height and tilt adjustment with backrest pivot hinge to provide free float tilting for lumbar comfort.