

SURGEONS FOOT OPERATED TILT CONTOUR STOOLS



CONTOUR SEATS

From **£390.00** Exc. VAT



HAN9640

HAN9645

Surgeons Foot Operated Tilt Contour Stools at a glance:

- Pneumatic height adjustment – standard 20"–28"
- Contoured shaped seat pad to aid correct posture
- Pneumatic 35° seat tilt adjustment
- 360° swivel action
- Anti-microbial vinyl covering options
- Hard free-wheeling castor wheels
- Available with Height adjustable torso/arm

Surgeons Foot Operated Tilt Contour Stool dimensions (standard):

Low and High height versions are available.

- Height: 22" – 27" (Adjustable)
- Wheel Base: 25"
- Max. User Weight 165kg

Production lead time: approx. 7-10 working days



UK Made. CE marked



Anti-microbial vinyl



Assembly required



FREE Delivery anywhere in mainland UK by courier

WITHOUT TORSO/ARM SUPPORT

HAN9640 | **£390.00** Exc. VAT

WITH TORSO/ARM SUPPORT

HAN9645 | **£490.00** Exc. VAT

Our Surgeon's Foot Operated Tilt Contour Stools take workplace comfort to a whole new level! These innovative contoured stools feature a shaped seat pad with 35° pneumatic tilt adjustment for further fine-tuning and posture management.

These Surgeons Foot Operated Tilt Contour Stools benefit from a 22" – 27" inch pneumatic height adjustment, hard, free-wheeling castor wheels and a choice of hard-wearing antimicrobial vinyl colour options.

This product is available both with and without an optional torso arm. The optional armrest is an ideal add-on for dentists undertaking regular time-consuming procedures and intricate, precision dental work.

Extra options:



Additional Memory Foam padding gives extra comfort if seated for long periods

+ £30



Wheel Variants:
Soft Anti-Static -roll freely on hard floors; Brake Loaded -roll & lock when sat on; Static -cup feet

+ £25



Height Variants:
Low Lift 18"-22"
High Lift 23"-31"

+ £20

Order online at meditelle-dental.co.uk or email sales@meditelleltd.co.uk

SURGEONS FOOT OPERATED TILT CONTOUR STOOLS

Specifications




HAN9640

HAN9645

UK Made. CE marked
 Anti-microbial vinyl
 Assembly required

FREE Delivery anywhere in mainland UK by courier

Technical information:



Height Information
 Standard 22"-27" (56cm - 69cm)
 Low option 18"-22" (46cm - 51cm)
 High option 23"-31" (58cm - 79cm)

SWL 165kg

Vinyl colours:



Special features:



Pneumatic Foot Operated Height Control – allows hands free, sterile height adjustment of the seat by foot.



Torso Arm Rest – Height adjustable swing support.
(HAN9645)