

# TILT CONTOUR CHAIRS

From **£340.00** Exc. VAT



HAN9650

HAN9655

## Tilt Contour Chairs at a glance:

- Pneumatic height adjustment – standard 20"–28"
- Contoured shaped seat pad to aid correct posture
- Contoured back support with dual adjustment
- Pneumatic 35° seat tilt adjustment
- 360° swivel action
- Anti-microbial vinyl covering options
- Hard free-wheeling castor wheels
- Available with and without a footrest ring

## Tilt Contour Chair dimensions (standard):

*Low and High height versions are available.*

- Height: 20" – 28" (Adjustable)
- Wheel Base: 25"
- Max. User Weight 165kg

**Production lead time: approx. 7-10 working days**



UK Made. CE marked



Anti-microbial vinyl



Assembly required



**FREE Delivery** anywhere in mainland UK by courier

**WITHOUT FOOT REST**

HAN9650 | **£340.00** Exc. VAT

**WITH FOOT REST**

HAN9655 | **£370.00** Exc. VAT

Our Tilt Contour Chairs are the height of workplace comfort! These innovative contoured chairs feature a shaped seat pad with 35° pneumatic tilt adjustment and a dual adjustment contoured backrest for further fine-tuning and posture management. These Tilt Contour Chairs benefit from a 20" – 28" inch pneumatic height adjustment, hard, free-wheeling castor wheels and a choice of hard-wearing antimicrobial vinyl colour options.

This dental chair is available both with and without an optional footrest. The optional footrest is an ideal add-on for dentists using the chair at full extension as the footrest ring helps to provide additional stability and comfort.

## Extra options:



Additional Memory Foam padding gives extra comfort if seated for long periods

**+ £30**



Wheel Variants:  
Soft Anti-Static –roll freely on hard floors; Brake Loaded –roll & lock when sat on; Static –cup feet

**+ £25**



Height Variants:  
Low Lift 19"–24"  
High Lift 22"–32"

**+ £20**

Order online at [meditelle-dental.co.uk](https://meditelle-dental.co.uk) or email [sales@meditelleltd.co.uk](mailto:sales@meditelleltd.co.uk)

# TILT CONTOUR CHAIRS

## Specifications



HAN9650




HAN9655

UK Made. CE marked
 Anti-microbial vinyl
 Assembly required

**FREE Delivery** anywhere in mainland UK by courier

### Technical information:



**Height Information**  
 Standard 20"-28" (51cm - 71cm)  
 Low option 19"-24" (48cm - 61cm)  
 High option 22"-32" (56cm - 81cm)

SWL 165kg

### Vinyl colours:



### Special features:



**Foot Rest Ring** - provides a foot rest support for added comfort and if working at a higher height.

(HAN9655)



**Contoured moulded back support** - features dual height and tilt adjustment with backrest pivot hinge to provide free float tilting for lumbar comfort.