

SURGEONS FOOT OPERATED TILT CONTOUR CHAIRS

From **£470.00** Exc. VAT



HAN9670

HAN9675

Surgeons Foot Operated Tilt Contour Chairs at a glance:

- Pneumatic height adjustment – standard 22”-27”
- Contoured shaped seat pad to aid correct posture
- Contoured back support with dual adjustment
- Pneumatic 35° seat tilt adjustment
- 360° swivel action
- Anti-microbial vinyl covering options
- Hard free-wheeling castor wheels
- Available with Height adjustable torso/arm

Surgeons Foot Operated Tilt Contour Chair dimensions (standard):

Low and High height versions are available.

- Height: 22” – 27” (Adjustable)
- Wheel Base: 25”
- Max. User Weight 165kg

Production lead time: approx. 7-10 working days

UK Made. CE marked Anti-microbial vinyl Assembly required

FREE Delivery anywhere in mainland UK by courier

WITHOUT TORSO/ARM SUPPORT HAN9670 | **£470.00** Exc. VAT

WITH TORSO/ARM SUPPORT HAN9675 | **£570.00** Exc. VAT

Our Surgeon's Foot Operated Tilt Contour Chairs offer comfort and stability when you need it most! These premium dentist's chairs feature an adjustable contoured seat pad and backrest, 35° seat tilt for posture management and a pneumatic foot-operated height adjustment which covers 22”-27” inches on our standard height model.

These Surgeon's Foot Operated Tilt Contour Chairs benefit from a full 360° range of motion, hard, free-wheeling castor wheels and a choice of antimicrobial vinyl upholstery options.

This dental chair is available with or without an option torso/armrest which can be useful for time-consuming or intricate dental procedures by providing additional stability and support to the surgeon's arm or torso.

Extra options:



Additional Memory Foam padding gives extra comfort if seated for long periods

+ £30



Wheel Variants:
Soft Anti-Static -roll freely on hard floors; Brake Loaded -roll & lock when sat on; Static -cup feet

+ £25



Height Variants:
Low Lift 18”-22”
High Lift 23”-31”

+ £20

Order online at meditelle-dental.co.uk or email sales@meditelleltd.co.uk

SURGEONS FOOT OPERATED TILT CONTOUR CHAIRS

Specifications




HAN9670

HAN9675

UK Made. CE marked Anti-microbial vinyl Assembly required

FREE Delivery anywhere in mainland UK by courier

Technical information:



Height Information
Standard 22"-27" (56cm - 69cm)
Low option 18"-22" (46cm - 51cm)
High option 23"-31" (58cm - 79cm)

SWL 165kg

Vinyl colours:



Special features:



Pneumatic Foot Operated Height Control - allows hands free, sterile height adjustment of the seat by foot.



Torso Arm Rest - Height adjustable swing support. (HAN9675)



Contoured moulded back support - features dual height and tilt adjustment with backrest pivot hinge to provide free float tilting for lumbar comfort.

Order online at meditelle-dental.co.uk or email sales@meditelleltd.co.uk