

# TILT CONTOUR CHAIRS WITH ARMS

From **£370.00** Exc. VAT



HAN9680

HAN9685

## Tilt Contour Chairs with Arms at a glance:

- Pneumatic height adjustment – standard 20”-28”
- Contoured shaped seat pad to aid correct posture
- Contoured back support with dual adjustment
- Pneumatic 35° seat tilt adjustment
- 360° swivel action
- Anti-microbial vinyl covering options
- Hard free-wheeling castor wheels
- Multi-adjusting side arm
- Available with and without a footrest ring

## Tilt Contour Chair with Arms dimensions (standard):

Low and High height versions are available.

- Height: 20” – 28” (Adjustable)
- Wheel Base: 25”
- Max. User Weight 165kg

**Production lead time: approx. 7-10 working days**

UK Made. CE marked Anti-microbial vinyl Assembly required

**FREE Delivery** anywhere in mainland UK by courier

**WITHOUT FOOT REST**

HAN9680 | **£370.00** Exc. VAT

**WITH FOOT REST**

HAN9685 | **£400.00** Exc. VAT

These stylish tilt contour chairs with arms are the ideal ergonomic seating solution for your dentist's office or surgery. Our tilt contour chairs with arms feature an adjustable contoured seat pad and backrest, 35° seat tilt adjustment and a pneumatic height adjustment that covers 20” – 28” inches on our standard height model.

The tilt contour chair with arms boasts a full 360° range of motion, dual multi-adjusting side arms, hard, free-wheeling castor wheels and a choice of anti-microbial vinyl upholstery options.

This innovative dentist's seat is available with or without an optional footrest ring. This useful add-on is perfectly suited to dentists who require additional foot support when working at the chair's maximum height elevation.

## Extra options:



Additional Memory Foam padding gives extra comfort if seated for long periods

**+ £30**



Wheel Variants:  
Soft Anti-Static –roll freely on hard floors; Brake Loaded –roll & lock when sat on; Static –cup feet

**+ £25**



Height Variants:  
Low Lift 19”-24”  
High Lift 22”-32”

**+ £20**

Order online at [meditelle-dental.co.uk](https://meditelle-dental.co.uk) or email [sales@meditelleltd.co.uk](mailto:sales@meditelleltd.co.uk)

# TILT CONTOUR CHAIRS WITH ARMS

## Specifications



HAN9680




HAN9685

UK Made. CE marked Anti-microbial vinyl Assembly required

**FREE Delivery** anywhere in mainland UK by courier

### Technical information:



**Height Information**  
 Standard 20"-28" (51cm - 71cm)  
 Low option 19"-24" (48cm - 61cm)  
 High option 22"-32" (56cm - 81cm)

SWL 165kg

### Vinyl colours:



### Special features:



**Foot Rest Ring** - provides a foot rest support for added comfort and if working at a higher height.

(HAN9685)



**Arm Rests** - Provide a multi adjustable arm support.



**Contoured moulded back support** - features dual height and tilt adjustment with backrest pivot hinge to provide free float tilting for lumbar comfort.

Order online at [meditelle-dental.co.uk](http://meditelle-dental.co.uk) or email [sales@meditelleltd.co.uk](mailto:sales@meditelleltd.co.uk)