

SURGEONS FOOT OPERATED TILT CONTOUR CHAIR WITH ARMS



CONTOUR SEATS

From **£500.00** Exc. VAT



HAN9690

Surgeons Foot Operated Tilt Contour Chair with Arms at a glance:

- Pneumatic height adjustment – standard 20"–28"
- Contoured shaped seat pad to aid correct posture
- Pneumatic 35° seat tilt adjustment
- 360° swivel action
- Anti-microbial vinyl covering options
- Hard free-wheeling castor wheels
- Multi-adjusting side arms

Surgeons Foot Operated Tilt Contour Chair with Arms dimensions (standard):

Low and High height versions are available.

- Height: 22" – 27" (Adjustable)
- Wheel Base: 25"
- Max. User Weight 165kg

Production lead time: approx. 7-10 working days



UK Made. CE marked



Anti-microbial vinyl



Assembly required



FREE Delivery anywhere in mainland UK by courier

HAN9690 | **£500.00** Exc. VAT

Our Surgeon's Foot Operated Tilt Contour Chairs with Arms take workplace comfort to a whole new level! These innovative contoured chairs feature a shaped seat pad with 35° pneumatic tilt adjustment for further fine-tuning and posture management and a comfortable dual adjustment backrest for added support.

These Surgeons Foot Operated Tilt Contour Chairs with Arms benefit from a 22" – 27" inch pneumatic height adjustment, hard, free-wheeling castor wheels, flexible, multi-adjusting side arms and a choice of hard-wearing antimicrobial vinyl colour options.

Extra options:



Additional Memory Foam padding gives extra comfort if seated for long periods

+ £30



Wheel Variants:
Soft Anti-Static -roll freely on hard floors; Brake Loaded -roll & lock when sat on; Static -cup feet

+ £25



Height Variants:
Low Lift 18"–22"
High Lift 23"–31"

+ £20

Order online at meditelle-dental.co.uk or email sales@meditelleltd.co.uk

SURGEONS FOOT OPERATED TILT CONTOUR CHAIR WITH ARMS

Specifications



HAN9690



UK Made. CE marked



Anti-microbial vinyl




Assembly required



FREE Delivery anywhere in mainland UK by courier

Technical information:



Height Information
Standard 22"-27" (56cm - 69cm)
Low option 18"-22" (46cm - 51cm)
High option 23"-31" (58cm - 79cm)



SWL 165kg

Vinyl colours:



Special features:



Pneumatic Foot Operated Height Control – allows hands free, sterile height adjustment of the seat by foot.



Arm Rests – Provide a multi adjustable arm support.



Contoured moulded back support – features dual height and tilt adjustment with backrest pivot hinge to provide free float tilting for lumbar comfort.