

DENTAL OPERATOR TUB STOOLS



HAN9700

HAN9710

Dental Tub Stools at a glance:

- Pneumatic height adjustment – standard 22”-30”
- 360° swivel action
- Anti-microbial vinyl covering
- Hard free-wheeling castor wheels
- Extra thick seat padding for maximum comfort
- Available with and without a footrest ring

Dental Tub Stool dimensions (standard):

Low and High height versions are available.

- Height: 22” – 30” (Adjustable)
- Seat Pad Diameter: 15.5”
- Wheel Base: 23”
- Max. User Weight 165kg

Production lead time: approx. 7-10 working days

 UK Made. CE marked  Anti-microbial vinyl  Assembly required

 **FREE Delivery** anywhere in mainland UK by courier

WITHOUT FOOT REST

HAN9700 | **£180.00** Exc. VAT

WITH FOOT REST

HAN9710 | **£210.00** Exc. VAT

Introducing our innovative Dental Tub Stools. These simple Tub Stools are the perfect addition to any dentist’s office or surgery. Our Dental Tub Stools are available in a range of anti-microbial vinyl colours and boast a practical, full 360° swivel design. Meditelle Dental Tub stools feature a smooth. Pneumatic height adjustment, hard, free-wheeling castor wheels and an extra thick seat pad for maximum comfort and durability.

This versatile tub stool is available with and without an additional footrest ring.

Extra options:



Additional Memory Foam padding gives extra comfort if seated for long periods

+ £30



Wheel Variants:
Soft Anti-Static –roll freely on hard floors; Brake Loaded –roll & lock when sat on; Static –cup feet

+ £25



Height Variants:
Low Lift 21”-26”
High Lift 24”-34”

+ £20

Order online at meditelle-dental.co.uk or email sales@meditelleltd.co.uk

DENTAL OPERATOR TUB STOOLS

Specifications




HAN9700

HAN9710

UK Made. CE marked Anti-microbial vinyl Assembly required

FREE Delivery anywhere in mainland UK by courier

Technical information:



Height Information
Standard 22"-30" (56cm - 76cm)
Low option 21"-26" (53cm - 66cm)
High option 24"-34" (61cm - 86cm)

SWL 165kg

Vinyl colours:



Special features:



Foot Rest Ring - provides a foot rest support for added comfort and if working at a higher height.

(HAN9710)

Order online at meditelle-dental.co.uk or email sales@meditelleltd.co.uk