

# DENTAL TUB CHAIRS WITH TILT

From **£295.00** Exc. VAT



HAN9720



HAN9730

## Dental Tub Chairs with Tilt at a glance:

- Pneumatic height adjustment – standard 22”-30”
- Pneumatic 35° seat tilt adjustment
- Contoured back support with dual adjustment
- 360° swivel action
- Anti-microbial vinyl covering
- Hard free-wheeling castor wheels
- Extra thick seat padding for maximum comfort
- Available with and without a footrest ring

## Dental Tub Chair with Tilt dimensions (standard):

Low and High height versions are available.

- Height: 22” – 30” (Adjustable)
- Seat Pad Diameter: 15.5”
- Wheel Base: 23”
- Max. User Weight 165kg

**Production lead time: approx. 7-10 working days**

 UK Made. CE marked  Anti-microbial vinyl  Assembly required

 **FREE Delivery** anywhere in mainland UK by courier

WITHOUT FOOT REST

HAN9720 | **£295.00** Exc. VAT

WITH FOOT REST

HAN9730 | **£325.00** Exc. VAT

Introducing our innovative Dental Tub Chairs with tilt and backrest. These creative Tub Chairs are the perfect addition to any dentist’s office or surgery. These seats boast a Contoured backrest with a comfortable Pneumatic 35° seat tilt adjustment, a pneumatic height adjustment and a full motion 360° swivel action.

Our Dental Tub Chairs with tilt are available in a choice of antimicrobial vinyl colours and feature extra thick seat padding to provide maximum comfort for long and intensive dental treatments.

The practical and comfortable dental stools are available with and without a footrest ring and both models benefit from hard, free-wheeling castor wheels.

## Extra options:



Additional Memory Foam padding gives extra comfort if seated for long periods

**+ £30**



Wheel Variants:  
Soft Anti-Static –roll freely on hard floors; Brake Loaded –roll & lock when sat on; Static –cup feet

**+ £25**



Height Variants:  
Low Lift 21”-26”  
High Lift 24”-34”

**+ £20**

Order online at [meditelle-dental.co.uk](https://meditelle-dental.co.uk) or email [sales@meditelleltd.co.uk](mailto:sales@meditelleltd.co.uk)

# DENTAL TUB CHAIRS WITH TILT

## Specifications



HAN9720



HAN9730

UK Made. CE marked
 Anti-microbial vinyl
 Assembly required

**FREE Delivery** anywhere in mainland UK by courier

### Technical information:



**Height Information**  
 Standard 22"-30" (56cm - 76cm)  
 Low option 21"-26" (53cm - 66cm)  
 High option 24"-34" (61cm - 86cm)

SWL 165kg

### Vinyl colours:



### Special features:



**Foot Rest Ring** - provides a foot rest support for added comfort and if working at a higher height.

(HAN9730)



**Contoured moulded back support** - features dual height and tilt adjustment with backrest pivot hinge to provide free float tilting for lumbar comfort.

Order online at [meditelle-dental.co.uk](http://meditelle-dental.co.uk) or email [sales@meditelleltd.co.uk](mailto:sales@meditelleltd.co.uk)