

DENTAL TUB STOOLS WITH TORSO ARM

From **£260.00** Exc. VAT



HAN9740

HAN9750

Dental Tub Stools with Torso Arm at a glance:

- Pneumatic height adjustment – standard 22”-30”
- Height adjustable torso/arm
- 360° swivel action
- Anti-microbial vinyl covering
- Hard free-wheeling castor wheels
- Extra thick seat padding for maximum comfort
- Available with and without a footrest ring

Dental Tub Stool with Torso Arm dimensions (standard):

Low and High height versions are available.

- Height: 22” – 30” (Adjustable)
- Seat Pad Diameter: 15.5”
- Wheel Base: 23”
- Max. User Weight 165kg

Production lead time: approx. 7-10 working days

 UK Made. CE marked  Anti-microbial vinyl  Assembly required

 **FREE Delivery** anywhere in mainland UK by courier

WITHOUT FOOT REST

HAN9740 | **£260.00** Exc. VAT

WITH FOOT REST

HAN9750 | **£290.00** Exc. VAT

Introducing our innovative Dental Tub Stools with Torso Arm. These creative Tub Stools are the perfect addition to any dentist’s office or surgery and the addition of a manoeuvrable torso arm allows the stool to be used for precision dental work and examinations.

Our Dental Tub Stools are available in a range of anti-microbial vinyl colours and boast a practical, full 360° swivel design. Meditelle Dental Tub stools feature a smooth. Pneumatic height adjustment, hard, free-wheeling castor wheels and an extra thick seat pad for maximum comfort and durability.

This versatile tub stool with Torso Arm is available with and without an additional footrest ring.

Extra options:



Additional Memory Foam padding gives extra comfort if seated for long periods

+ £30



Wheel Variants:
Soft Anti-Static –roll freely on hard floors; Brake Loaded –roll & lock when sat on; Static –cup feet

+ £25



Height Variants:
Low Lift 21”-26”
High Lift 24”-34”

+ £20

Order online at meditelle-dental.co.uk or email sales@meditelleltd.co.uk

DENTAL TUB STOOLS WITH TORSO ARM

Specifications



HAN9740




HAN9750

UK Made. CE marked Anti-microbial vinyl Assembly required

FREE Delivery anywhere in mainland UK by courier

Technical information:



Height Information
 Standard 22"-30" (56cm - 76cm)
 Low option 21"-26" (53cm - 66cm)
 High option 24"-34" (61cm - 86cm)

SWL 165kg

Vinyl colours:



Special features:



Foot Rest Ring - provides a foot rest support for added comfort and if working at a higher height.

(HAN9750)



Torso Arm Rest - Height adjustable swing support.