

SHAPED TUB CHAIRS WITH TILT

From **£330.00** Exc. VAT



HAN9830

HAN9840

Shaped Tub Chairs with Tilt at a glance:

- Pneumatic height adjustment – standard 22”-30”
- Pneumatic 35° seat tilt adjustment
- Contoured back support with dual adjustment
- 360° swivel action
- Anti-microbial vinyl covering
- Hard free-wheeling castor wheels
- Extra thick seat padding for maximum comfort
- Available with and without a footrest ring

Shaped Tub Chair with Tilt dimensions (standard):

Low and High height versions are available.

- Height: 22” – 30” (Adjustable)
- Seat Pad Diameter: 15.5”
- Wheel Base: 23”
- Max. User Weight 165kg

Production lead time: approx. 7-10 working days

UK Made. CE marked Anti-microbial vinyl Assembly required

FREE Delivery anywhere in mainland UK by courier

WITHOUT FOOT REST

HAN9830 | **£330.00** Exc. VAT

WITH FOOT REST

HAN9840 | **£360.00** Exc. VAT

Introducing our innovative Shaped Tub Chairs with tilt and backrest. These creative Tub Chairs are the perfect addition to any dentist's office or surgery. These seats boast a Contoured backrest with a comfortable Pneumatic 35° seat tilt adjustment, a pneumatic height adjustment and a full motion 360° swivel action.

Our Shaped Tub Chairs with tilt are available in a choice of antimicrobial vinyl colours and feature extra thick seat padding to provide maximum comfort for long and intensive dental treatments.

These practical and comfortable shaped chairs are available with and without a footrest ring and both models benefit from hard, free-wheeling castor wheels.

Extra options:



Additional Memory Foam padding gives extra comfort if seated for long periods

+ £30



Wheel Variants:
Soft Anti-Static –roll freely on hard floors; Brake Loaded –roll & lock when sat on; Static –cup feet

+ £25



Height Variants:
Low Lift 21”-26”
High Lift 24”-34”

+ £20

Order online at [meditelle-dental.co.uk](https://www.meditelle-dental.co.uk) or email sales@meditelleltd.co.uk

SHAPED TUB CHAIRS WITH TILT

Specifications



HAN9830




HAN9840

UK Made. CE marked Anti-microbial vinyl Assembly required

FREE Delivery anywhere in mainland UK by courier

Technical information:



Height Information
 Standard 22"-30" (56cm - 76cm)
 Low option 21"-26" (53cm - 66cm)
 High option 24"-34" (61cm - 86cm)

SWL 165kg

Vinyl colours:



Special features:



Foot Rest Ring - provides a foot rest support for added comfort and if working at a higher height.

(HAN9840)



Contoured moulded back support - features dual height and tilt adjustment with backrest pivot hinge to provide free float tilting for lumbar comfort.

Order online at meditelle-dental.co.uk or email sales@meditelleltd.co.uk