

# SHAPED TUB STOOLS WITH TORSO ARM

From **£295.00** Exc. VAT



HAN9860

HAN9870

## Shaped Tub Stools with Torso Arm at a glance:

- Pneumatic height adjustment – standard 22”-30”
- Operator Stools with extra thick seat pad
- Shaped seat prevents upper leg impingement
- 360° swivel action
- Anti-microbial vinyl covering
- Hard free-wheeling castor wheels
- Height adjustable torso/arm
- Extra thick seat padding for maximum comfort
- Available with and without a footrest ring

## Shaped Tub Stool with Torso Arm dimensions (standard):

*Low and High height versions are available.*

- Height: 22” – 30” (Adjustable)
- Seat Pad Diameter: 15.5”
- Wheel Base: 23”
- Max. User Weight 165kg

**Production lead time: approx. 7-10 working days**

UK Made. CE marked Anti-microbial vinyl Assembly required

**FREE Delivery** anywhere in mainland UK by courier

**WITHOUT FOOT REST**

HAN9860 | **£295.00** Exc. VAT

**WITH FOOT REST**

HAN9870 | **£325.00** Exc. VAT

Our Shaped Tub Stools with Torso Arms are built for maximum comfort and stability making them the ideal seating solution for any dentist's office or surgery. These functional tub stools boast an extra thick seat pad as standard, and a shaped seat to prevent upper leg impingement while hard, free-wheeling castors keep you firmly rooted to the floor. This shaped tool stool is available in a wide range of antimicrobial vinyl colour options and boasts a pneumatic height adjustment and 360° swivel action as well as a manoeuvrable torso arm for additional support when performing intricate dental work.

*Our Meditelle Shaped Tub Stools with Torso Arms are available with and without a footrest ring.*

## Extra options:



Additional Memory Foam padding gives extra comfort if seated for long periods

**+ £30**



Wheel Variants:  
Soft Anti-Static –roll freely on hard floors; Brake Loaded –roll & lock when sat on; Static –cup feet

**+ £25**



Height Variants:  
Low Lift 21”-26”  
High Lift 24”-34”

**+ £20**

Order online at [meditelle-dental.co.uk](https://meditelle-dental.co.uk) or email [sales@meditelleltd.co.uk](mailto:sales@meditelleltd.co.uk)

# SHAPED TUB STOOLS WITH TORSO ARM

## Specifications




HAN9860

HAN9870

UK Made. CE marked   Anti-microbial vinyl   Assembly required

**FREE Delivery** anywhere in mainland UK by courier

### Technical information:



**Height Information**  
 Standard 22"-30" (56cm - 76cm)  
 Low option 21"-26" (53cm - 66cm)  
 High option 24"-34" (61cm - 86cm)

SWL 165kg

### Vinyl colours:



### Special features:



**Foot Rest Ring** - provides a foot rest support for added comfort and if working at a higher height.

(HAN9870)



**Torso Arm Rest** - Height adjustable swing support.

Order online at [meditelle-dental.co.uk](http://meditelle-dental.co.uk) or email [sales@meditelleltd.co.uk](mailto:sales@meditelleltd.co.uk)