

## Antimicrobial Vinyl Information

Meditelle Vinyl Fabrics include Biocides, an antimicrobial ingredient which ensures it is safe for use in a medical environment.

At Meditelle we use a vinyl that contains Sanitized® antimicrobial formula which has been specifically developed to provide a high level of immediate and long-term protection against attack by micro-organisms. The Sanitized® hygiene function reduces the growth of micro-organisms such as Clostridia and antibiotic resistant bacteria such as MRSA, VRE and NDM-1.

**The Biocide ingredient locked within our fabric also deactivates and prevents the growth of Coronavirus, it will eliminate the virus in just over 1 hour of contact without any extra disinfecting or cleaning. Tested to ISO 18184:2019.**

The Sanitized® technology is added directly to the high-performance vinyl we use during manufacture, protecting from bacteria around the clock. This is especially important on surfaces that are frequently contaminated with germs and where ongoing protection between cleaning cycles is indispensable.

No matter how many times our vinyl is washed or cleaned the antimicrobial ingredients will remain active for the life of the fabric. The inherent antimicrobial properties of the vinyls we use ensures they have an unparalleled safety profile.

Sanitized® will not leach into the environment, nor transfer to other articles or to the skin. It prevents bacterial and fungal growth as well as eliminating microbial staining.



Independent UK laboratory tested.

Test Certificate copies available upon request.

**Meditelle Ltd**

**T:** 0121 332 1850

**E:** [enquiry@meditelle.co.uk](mailto:enquiry@meditelle.co.uk)

**W:** [meditelle.co.uk](http://meditelle.co.uk) • [meditelle-dental.co.uk](http://meditelle-dental.co.uk)